



Product Specification
Micro sample tube, 1.000 µl haemolysis solution/20 µl, 20 µl, push cap

Page 1



Product description

Order number	41.1438.904
Product description	Micro sample tube, 1.000 µl haemolysis solution/20 µl, capillary blood collection, preparation: SARSTEDT haemolysis solution, incl. 20 µl end-to-end capillary, nominal volume: 20 µl, (LxØ) with cap: 42 x 10.8 mm, cap: yellow, label/print: black, push cap, conical base, 1, 100 piece(s)/bag

Product characteristics

Type of collection	capillary
Type of preparation	SARSTEDT haemolysis solution
Label/ Print	with print
Colour of print/label	black
Cap	flat cap
Closure type	push cap
Application	capillary blood collection
Example of use	blood sugar determination from capillary blood
Base shape	conical base
Volume of work	20 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Micro sample tube, 1.000 µl haemolysis
solution/20 µl, 20 µl, push cap**

Page 2

Size

Sample volume	20 µl
Diameter	10.8 mm
Length including cap	42 mm
Length excluding cap	39 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	yellow

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	96000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0025 kg
Weight of case	2.81 kg
EAN of inner box	4038917129169
EAN case	4038917022002

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

