



**Product Specification**  
**Micro sample tube, 1,000 µl ESAT/20µl, 20 µl,  
push cap**

**Page 1**



**Product description**

Order number	41.1443.904
Product description	Micro sample tube, 1,000 µl ESAT/20µl, capillary blood collection, preparation: ESAT haemolysis solution, incl. 20 µl end-to-end capillary, nominal volume: 20 µl, (LxØ) with cap: 42 x 10.8 mm, cap: yellow, push cap, skirted conical base, 2, 100 piece(s)/bag

**Product characteristics**

Type of collection	capillary
Type of preparation	ESAT haemolysis solution
Cap	flat cap
Closure type	push cap
Application	capillary blood collection
Example of use	Blood sugar determination from capillary blood
Base shape	skirted conical base
Volume of work	20 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro sample tube, 1,000 µl ESAT/20µl, 20 µl,  
push cap**

**Page 2**

**Size**

Sample volume	20 µl
Diameter	10.8 mm
Length including cap	42 mm
Length excluding cap	39 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	yellow

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	96000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0026 kg
Weight of case	2.93 kg
EAN of inner box	4038917129237
EAN case	4038917075985

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**