



Product description

Order number	42.117
Product description	Micro sample tube, preparation: 1% brilliant cresyl blue solution, nominal volume: 100 µl, (LxØ) with cap: 46 x 10.8 mm, cap: blue, screw cap, conical base

Product characteristics

Type of collection	Venous and capillary
Type of preparation	1% brilliant cresyl blue solution
Preparation concentration	1% brilliant cresyl blue solution
Type of cap	screw cap
Example of use	Reticulocyte counting
Base shape	conical base
Volume of work	100 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Micro sample tube, Brilliant cresyl blue, 100 µl,
screw cap**

Page 2

Size

Sample volume	100 µl
Diameter	10,8 mm
Length of product	46 mm
Length including cap	46 mm
Length excluding cap	44 mm

Material & colours

Colour of product	transparent
Cap material	Polypropylene (PP)
Colour of cap	blue

Purity & certification

Product category	IVD product
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	inner box
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	117000
Depth of case	283 mm
Width of case	173 mm
Height of case	135 mm
Case volume	0,0066 cbm
Weight of product	0,0015 kg
Weight of case	0,92 kg
EAN of inner box	4038917080897
EAN case	4038917365130

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-24

This document was prepared by EDP support and is valid without signature.



SARSTEDT