



**Product Specification**  
**Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent**

**Page 1**



**Product description**

Order number 80.620  
Product description Faeces tube, with spoon, transparent, push cap, cap: brown, (LxØ) without cap: 75 x 23.5 mm, 50 piece(s)/bag

**Product characteristics**

Closure type push cap  
Design spoon

**Size**

Diameter 23.5 mm  
Length of product 87 mm  
Length excluding cap 75 mm

**Material & colours**

Colour of product transparent  
Closure material Low Density Polyethylene (LD-PE)  
Colour of cap brown  
Tube material Polystyrene (PS)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent**

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	18000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0091 kg
Weight of case	5.29 kg
EAN of inner box	4038917082952
EAN case	4038917001939

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

