



**Product Specification**  
**Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent, sterile**

**Page 1**



**Product description**

Order number	80.621
Product description	Faeces tube, with spoon, transparent, push cap, cap: brown, (LxØ) without cap: 75 x 23.5 mm, sterile, 50 piece(s)/bag

**Product characteristics**

Closure type	push cap
Design	spoon

**Size**

Diameter	23.5 mm
Length of product	87 mm
Length excluding cap	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent, sterile**

**Page 2**

### Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	brown
Tube material	Polystyrene (PS)

### Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	13000
Depth of case	548 mm
Width of case	380 mm
Height of case	130 mm
Case volume	0.0271 cbm
Weight of product	0.0091 kg
Weight of case	2.73 kg
EAN of inner box	4038917082969
EAN case	4038917001946

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**