



**Product Specification**  
**Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent**

**Page 1**



**Product description**

|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 80.622.111                                                                                                                                          |
| Product description | Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon, enclosed, cap: brown, (LxØ) without cap: 107 x 25 mm, 250 piece(s)/bag |

**Product characteristics**

|              |                            |
|--------------|----------------------------|
| Cap          | Screw cap with stool spoon |
| Closure type | screw cap                  |
| Design       | spoon                      |
| Base shape   | flat-conical base          |

**Size**

|                      |        |
|----------------------|--------|
| Diameter             | 25 mm  |
| Length of product    | 111 mm |
| Length excluding cap | 107 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent**

**Page 2**

**Material & colours**

|                   |                                   |
|-------------------|-----------------------------------|
| Colour of product | transparent                       |
| Closure material  | High Density Polyethylene (HD-PE) |
| Colour of cap     | brown                             |
| Stem material     | Polystyrene (PS)                  |
| Tube material     | Polypropylene (PP)                |

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 250           |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 12000         |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 415 mm        |
| Case volume                   | 0.0657 cbm    |
| Weight of product             | 0.0129 kg     |
| Weight of case                | 7.14 kg       |
| EAN of inner box              | 4038917115346 |
| EAN case                      | 4038917014533 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**