



## Product Specification

Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent

Page 1



### Product description

Order number	80.623.022
Product description	Faeces tube, with spoon, transparent, screw cap, cap: brown, (LxØ) without cap: 101 x 16.5 mm, with paper label, 250 piece(s)/bag

### Product characteristics

Label/ Print	with paper label
Colour of print/label	black
Closure type	screw cap
Design	spoon

### Size

Diameter	16.5 mm
Length of product	107 mm
Length excluding cap	101 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent**

**Page 2**

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

### Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0079 kg
Weight of case	8.62 kg
EAN of inner box	4038917115407
EAN case	4038917014601

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**