




**Product Specification**  
**Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent**

**Page 1**



### Product description

|                     |   |
|---------------------|---|
| Order number        | 80.623.111  |
| Product description | Faeces tube, with spoon, transparent, screw cap, cap: brown, (LxØ) without cap: 101 x 16.5 mm, 500 piece(s)/bag |

### Product characteristics

|              |           |
|--------------|-----------|
| Closure type | screw cap |
| Design       | spoon     |

### Size

|                      |         |
|----------------------|---------|
| Diameter             | 16.5 mm |
| Length of product    | 107 mm  |
| Length excluding cap | 101 mm  |

### Material & colours

|                   |                                   |
|-------------------|-----------------------------------|
| Colour of product | transparent                       |
| Closure material  | High Density Polyethylene (HD-PE) |
| Colour of cap     | brown                             |
| Tube material     | Polypropylene (PP)                |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent**

**Page 2**

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 22000         |
| Depth of case                 | 568 mm        |
| Width of case                 | 348 mm        |
| Height of case                | 155 mm        |
| Case volume                   | 0.0306 cbm    |
| Weight of product             | 0.0068 kg     |
| Weight of case                | 3.89 kg       |
| EAN of inner box              | 4038917115520 |
| EAN case                      | 4038917014748 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-22

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**