




**Product Specification**  
**Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent**

**Page 1**



### Product description

Order number	80.623.111
Product description	Faeces tube, with spoon, transparent, screw cap, cap: brown, (LxØ) without cap: 101 x 16.5 mm, 500 piece(s)/bag

### Product characteristics

Closure type	screw cap
Design	spoon

### Size

Diameter	16.5 mm
Length of product	107 mm
Length excluding cap	101 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent**

**Page 2**

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	22000
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0.0306 cbm
Weight of product	0.0068 kg
Weight of case	3.89 kg
EAN of inner box	4038917115520
EAN case	4038917014748

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**