




Product Specification
Faeces tube, screw cap, (LxØ): 76 x 20 mm,
transparent, sterile

Page 1



Product description

Order number	80.734.311
Product description	Faeces tube, transparent, screw cap, cap: brown, (LxØ) without cap: 76 x 20 mm, with paper label, sterile

Product characteristics

Label/ Print	with paper label
Colour of print/label	black
Closure type	screw cap

Size

Diameter	20 mm
Length of product	77 mm
Length excluding cap	76 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-01-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Faeces tube, screw cap, (LxØ): 76 x 20 mm,
transparent, sterile

Page 2

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	20000
Depth of case	598 mm
Width of case	398 mm
Height of case	165 mm
Case volume	0.0393 cbm
Weight of product	0.0071 kg
Weight of case	4.49 kg
EAN of inner box	4038917119962
EAN case	4038917016964

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-01-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT