



Product Specification
S-Monovette® rack D17, Ø opening: 17 mm, 5 x 10, white



Product description

Order number	93.852.170
Product description	S-Monovette® rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, white, material: PP

Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	17 mm
Width of product	109 mm
Height of product	45 mm
Length of product	209 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
S-Monovette® rack D17, Ø opening: 17 mm, 5 x 10, white

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product white

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging -
Piece(s) / inner box 1
Piece (s) / outer case 20
Piece(s) / pallet 760
Depth of case 568 mm
Width of case 433 mm
Height of case 120 mm
Case volume 0.0295 cbm
Weight of product 0.1958 kg
Weight of case 4.62 kg
EAN of inner box 4038917122054
EAN case 4038917044905

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-11-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT