



Product description

Order number	93.852.171
Product description	Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, red, material: recycled PP

Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	17 mm
Width of product	109 mm
Height of product	45 mm
Length of product	209 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material recycled Polypropylene (recycled PP)
Colour of product red

Purity & certification

Product category not regulated
CE certified no CE mark
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging -
Piece(s) / inner box 1
Piece (s) / outer case 20
Piece(s) / pallet 760
Depth of case 568 mm
Width of case 433 mm
Height of case 120 mm
Case volume 0.0295 cbm
Weight of product 0.1958 kg
Weight of case 4.62 kg
EAN of inner box 4038917122061
EAN case 4038917018128

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

