



**Product Specification**  
**Block Rack D17, Ø opening: 17 mm, 5 x 10, orange**

**Page 1**



### Product description

Order number	93.852.175
Product description	Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, orange, material: PP

### Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

### Size

Diameter of opening	17 mm
Width of product	109 mm
Height of product	45 mm
Length of product	209 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Block Rack D17, Ø opening: 17 mm, 5 x 10,  
orange**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	orange

**Purity & certification**

Product category	not regulated
CE certified	no CE mark
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	760
Depth of case	568 mm
Width of case	433 mm
Height of case	120 mm
Case volume	0.0295 cbm
Weight of product	0.1958 kg
Weight of case	4.62 kg
EAN of inner box	4038917122108
EAN case	4038917018166

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**