




**Product Specification**  
**Block Rack D17, Ø opening: 17 mm, 5 x 10,  
orange**

**Page 1**



**Product description**

|                     |                                                                                                 |
|---------------------|-------------------------------------------------------------------------------------------------|
| Order number        | 93.852.175                                                                                      |
| Product description | Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, orange, material: PP |

**Product characteristics**

|                        |            |
|------------------------|------------|
| Label/ Print           | with print |
| Number of wells        | 50         |
| Storage capacity       | 50         |
| Format                 | 5 x 10     |
| Autoclavable (max. °C) | 121 °C     |

**Size**

|                     |        |
|---------------------|--------|
| Diameter of opening | 17 mm  |
| Width of product    | 109 mm |
| Height of product   | 45 mm  |
| Length of product   | 209 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Block Rack D17, Ø opening: 17 mm, 5 x 10,  
orange**

**Page 2**

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | orange             |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | no                      |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1             |
| Type of smallest subpackaging | -             |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 20            |
| Piece(s) / pallet             | 760           |
| Depth of case                 | 568 mm        |
| Width of case                 | 433 mm        |
| Height of case                | 120 mm        |
| Case volume                   | 0.0295 cbm    |
| Weight of product             | 0.1958 kg     |
| Weight of case                | 4.62 kg       |
| EAN of inner box              | 4038917122108 |
| EAN case                      | 4038917018166 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**