




**Product Specification**  
**Double block rack D17, Ø opening: 17 mm, 10 x 10, grey**

**Page 1**



**Product description**

Order number	93.852.200
Product description	Double block rack D17, for 100 collection tubes, Ø opening: 17 mm, format: 10 x 10, grey, material: PP

**Product characteristics**

Number of wells	100
Storage capacity	100
Format	10 x 10
Autoclavable (max. °C)	121 °C

**Size**

Diameter of opening	17 mm
Width of product	218 mm
Height of product	45 mm
Length of product	209 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Double block rack D17, Ø opening: 17 mm, 10 x 10, grey**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	grey

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	12
Piece(s) / pallet	528
Depth of case	453 mm
Width of case	433 mm
Height of case	165 mm
Case volume	0.0324 cbm
Weight of product	0.3917 kg
Weight of case	5.41 kg
EAN of inner box	4038917122146
EAN case	4038917243124

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

