



### **Product description**

Order number	93.852.300
Product description	Double block rack D17, for 100 collection tubes, Ø opening: 17 mm, format: 5 x 20, grey, material: recycled PP

### **Product characteristics**

Number of wells	100
Storage capacity	100
Format	5 x 20
Autoclavable (max. °C)	121 °C

### **Size**

Diameter of opening	17 mm
Width of product	109 mm
Height of product	45 mm
Length of product	418 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material recycled Polypropylene (recycled PP)  
Colour of product grey

### **Purity & certification**

Product category not regulated  
CE certified kein CE-Kennzeichen  
Batched no

### **Packaging**

Minimum order qty. 1  
Type of smallest subpackaging -  
Piece(s) / inner box 1  
Piece (s) / outer case 12  
Piece(s) / pallet 528  
Depth of case 453 mm  
Width of case 433 mm  
Height of case 165 mm  
Case volume 0.0324 cbm  
Weight of product 0.3917 kg  
Weight of case 5.41 kg  
EAN of inner box 4038917122221  
EAN case 4038917018296

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

