




**Product Specification**  
**Double block rack D17, Ø opening: 17 mm, 10 x 10, grey, connected at head end**

**Page 1**



### Product description

Order number	93.852.300
Product description	Double block rack D17, for 100 collection tubes, Ø opening: 17 mm, format: 10 x 10, grey, connected at head end, material: PP

### Product characteristics

Number of wells	100
Storage capacity	100
Format	10 x 10
Autoclavable (max. °C)	121 °C

### Size

Diameter of opening	17 mm
Width of product	218 mm
Height of product	45 mm
Length of product	209 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-06-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Double block rack D17, Ø opening: 17 mm, 10 x 10, grey, connected at head end**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	grey

**Purity & certification**

Product in accordance with MDD or IVDD	no medical device   IVD
Lot or serial number	no

**Packaging**

Minimum order qty.	1
Type of inner packaging	inner box
Piece(s) / inner box	1
Piece(s) / case	12
Piece(s) / pallet	528
Depth of case	453 mm
Width of case	433 mm
Height of case	165 mm
Case volume	0,0324 cbm
Weight of product	0,3917 kg
Weight of case	5,41 kg
EAN of inner box	4038917122221
EAN case	4038917018296

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

