



Product description

Order number	93.853
Product description	Block Rack D13, for 50 Monovettes and tubes up to 13 mm diameter, Ø opening: 13 mm, format: 5 x 10, grey, material: recycled PP

Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	13 mm
Width of product	96 mm
Height of product	45 mm
Length of product	180 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	recycled Polypropylene (recycled PP)
Colour of product	grey

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	16
Piece(s) / pallet	1440
Depth of case	398 mm
Width of case	368 mm
Height of case	120 mm
Case volume	0.0176 cbm
Weight of product	0.1756 kg
Weight of case	3.27 kg
EAN of inner box	4038917084390
EAN case	4038917074018

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

