



Product Specification
S-Monovette® rack D13, Ø opening: 13 mm, 5 x 10, yellow



Product description

Order number	93.853.132
Product description	S-Monovette® rack D13, for 50 collection tubes, Ø opening: 13 mm, format: 5 x 10, yellow, material: PP

Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	13 mm
Width of product	96 mm
Height of product	45 mm
Length of product	180 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
S-Monovette® rack D13, Ø opening: 13 mm, 5 x 10, yellow

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product yellow

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of inner packaging -
Piece(s) / inner box 1
Piece(s) / case 16
Piece(s) / pallet 1440
Depth of case 398 mm
Width of case 368 mm
Height of case 120 mm
Case volume 0.0176 cbm
Weight of product 0.1687 kg
Weight of case 3.16 kg
EAN of inner box 4038917122535
EAN case 4038917074032

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-05-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT