



Product Specification
**Block Rack D13, Ø opening: 13 mm, 5 x 10,
orange**

Page 1



Product description

| | |
|---------------------|--|
| Order number | 93.853.135 |
| Product description | Block Rack D13, for 50 collection tubes, Ø opening: 13 mm, format: 5 x 10, orange, material: PP |

Product characteristics

| | |
|------------------------|------------|
| Label/ Print | with print |
| Number of wells | 50 |
| Storage capacity | 50 |
| Format | 5 x 10 |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|---------------------|--------|
| Diameter of opening | 13 mm |
| Width of product | 96 mm |
| Height of product | 45 mm |
| Length of product | 180 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Block Rack D13, Ø opening: 13 mm, 5 x 10,
orange**

Page 2

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | orange |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | - |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 16 |
| Piece(s) / pallet | 1440 |
| Depth of case | 398 mm |
| Width of case | 368 mm |
| Height of case | 120 mm |
| Case volume | 0.0176 cbm |
| Weight of product | 0.1687 kg |
| Weight of case | 3.16 kg |
| EAN of inner box | 4038917122566 |
| EAN case | 4038917074063 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT