



### **Product description**

Order number	93.854.122
Product description	S-Monovette® rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 10 x 5, yellow, material: PP

### **Product characteristics**

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	10 x 5
Autoclavable (max. °C)	121 °C

### **Size**

Diameter of opening	12 mm
Width of product	80 mm
Height of product	45 mm
Length of product	152 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow

### **Purity & certification**

Product in accordance with MDD or IVDD	no medical device   IVD
Lot or serial number	no

### **Packaging**

Minimum order qty.	1
Type of inner packaging	inner box
Piece(s) / inner box	1
Piece(s) / case	20
Piece(s) / pallet	1560
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm
Case volume	0,018 cbm
Weight of product	0,1214 kg
Weight of case	2,8 kg
EAN of inner box	4038917122658
EAN case	4038917018494

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

