



Product description

| | |
|---------------------|---|
| Order number | 93.854.122 |
| Product description | Block Rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 5 x 10, yellow, material: PP |

Product characteristics

| | |
|------------------------|------------|
| Label/ Print | with print |
| Number of wells | 50 |
| Storage capacity | 50 |
| Format | 5 x 10 |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|---------------------|--------|
| Diameter of opening | 12 mm |
| Width of product | 80 mm |
| Height of product | 45 mm |
| Length of product | 152 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product yellow

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging -
Piece(s) / inner box 1
Piece (s) / outer case 20
Piece(s) / pallet 1560
Depth of case 428 mm
Width of case 323 mm
Height of case 130 mm
Case volume 0.018 cbm
Weight of product 0.1214 kg
Weight of case 2.8 kg
EAN of inner box 4038917122658
EAN case 4038917018494

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

