



## Product description

Order number	93.854.124
Product description	Block Rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 5 x 10, blue, material: PP

## Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

## Size

Diameter of opening	12 mm
Width of product	80 mm
Height of product	45 mm
Length of product	152 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product blue

### **Purity & certification**

Product category not regulated  
CE certified no CE mark  
Batched no

### **Packaging**

Minimum order qty. 1  
Type of smallest subpackaging -  
Piece(s) / inner box 1  
Piece (s) / outer case 20  
Piece(s) / pallet 1560  
Depth of case 390 mm  
Width of case 355 mm  
Height of case 173 mm  
Case volume 0.024 cbm  
Weight of product 0.1214 kg  
Weight of case 2.8 kg  
EAN of inner box 4038917122672  
EAN case 4038917018517

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

