



Product Specification
Block Rack D12, Ø opening: 12 mm, 5 x 10, orange

Page 1



Product description

Order number	93.854.125
Product description	Block Rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 5 x 10, orange, material: PP

Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	12 mm
Width of product	80 mm
Height of product	45 mm
Length of product	152 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Block Rack D12, Ø opening: 12 mm, 5 x 10,
orange**

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product orange

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging -
Piece(s) / inner box 1
Piece (s) / outer case 20
Piece(s) / pallet 1560
Depth of case 428 mm
Width of case 323 mm
Height of case 130 mm
Case volume 0.018 cbm
Weight of product 0.1214 kg
Weight of case 2.8 kg
EAN of inner box 4038917122689
EAN case 4038917018524

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT