



Product description

Order number	93.854
Product description	S-Monovette® rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 10 x 5, grey, material: PP

Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	10 x 5
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	12 mm
Width of product	80 mm
Height of product	45 mm
Length of product	152 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
S-Monovette® rack D12, Ø opening: 12 mm, 10 x 5, grey

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product grey

Purity & certification

Product in accordance with MDD or IVDD no medical device | IVD
Lot or serial number no

Packaging

Minimum order qty. 1
Type of inner packaging inner box
Piece(s) / inner box 1
Piece(s) / case 20
Piece(s) / pallet 1560
Depth of case 428 mm
Width of case 323 mm
Height of case 130 mm
Case volume 0,018 cbm
Weight of product 0,1214 kg
Weight of case 2,8 kg
EAN of inner box 4038917084406
EAN case 4038917003247

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-04-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT