



## Product Specification

Tube with false bottom, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, flat tube bottom, PP, without  
cap, 100 piece(s)/bag

Page 1



### Product description

Order number	60.611.010
Product description	Tube with false bottom, working volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, flat tube bottom, transparent, material: PP, without cap, 100 piece(s)/bag, 1,000 piece(s)/carton

### Product characteristics

Graduation	no
Base shape	conical false bottom, flat tube bottom
Volume of work	5 ml

### Size

Diameter	15.3 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Tube with false bottom, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, flat tube bottom, PP, without  
cap, 100 piece(s)/bag**

**Page 2**

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	310 mm
Case volume	0.0491 cbm
Weight of product	0.004 kg
Weight of case	4.68 kg
EAN of inner box	4038917114820
EAN case	4038917014427

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-15

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**