



**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm, flat false bottom, PP, 100 piece(s)/bag**

**Page 1**



### Product description

Order number	60.611.310
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 92 x 15.3 mm, flat false bottom, transparent, material: PP, without cap, 100 piece(s)/bag

### Product characteristics

Base shape	flat false bottom
Volume of work	5 ml
Shape of inner tube	conical

### Size

Diameter	15,3 mm
Height of product	92 mm
Length of product	92 mm
Length excluding cap	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-04-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm, flat  
false bottom, PP, 100 piece(s)/bag**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Tube material	Polypropylene (PP)

**Purity & certification**

Product in accordance with MDD or IVDD	no medical device   IVD
Lot or serial number	yes

**Packaging**

Minimum order qty.	1000
Type of inner packaging	bag
Piece(s) / inner box	1000
Piece(s) / case	1000
Piece(s) / pallet	30000
Case volume	0,046 cbm
Weight of product	0,0041 kg
Weight of case	4,55 kg
EAN of inner box	4038917559799
EAN case	4038917367905

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-04-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**