



**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm, flat false bottom, PP, 100 piece(s)/bag**



**Product description**

Order number	62.611.300
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 92 x 15.3 mm, flat false bottom, transparent, material: PP, with print, print white, with graduations, cap assembled, 100 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Base shape	flat false bottom
Volume of work	5 ml
Shape of inner tube	conical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm, flat false bottom, PP, 100 piece(s)/bag**

**Size**

Diameter	15,3 mm
Height of product	92 mm
Length of product	92 mm
Length excluding cap	92 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Tube material	Polypropylene (PP)

**Purity & certification**

Product in accordance with MDD or IVDD	no medical device   IVD
Lot or serial number	yes

**Packaging**

Minimum order qty.	1000
Type of inner packaging	bag
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	36000
Case volume	0,0492 cbm
Weight of product	0,005 kg
Weight of case	5,56 kg
EAN of inner box	4038917368216
EAN case	4038917368209

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

