



## Product Specification

**Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm,  
conical false bottom, flat tube bottom, PP, cap  
assembled, 100 piece(s)/bag**

**Page 1**



### Product description

Order number	62.613
Product description	Screw cap tube, working volume: 3.5 ml, (LxØ): 92 x 13 mm, conical false bottom, flat tube bottom, transparent, material: PP, with graduations, cap assembled, 100 piece(s)/bag, 1,000 piece(s)/carton

### Product characteristics

Graduation	yes
Base shape	conical false bottom, flat tube bottom
Volume of work	3,5 ml

### Size

Diameter	13 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm,  
conical false bottom, flat tube bottom, PP, cap  
assembled, 100 piece(s)/bag**

**Page 2**

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	20000
Depth of box	340 mm
Width of box	330 mm
Height of box	20 mm
Depth of case	460 mm
Width of case	410 mm
Height of case	210 mm
Case volume	0.0396 cbm
Weight of product	0.0033 kg
Weight of case	4.34 kg
EAN of inner box	4038917082891
EAN case	4038917079136

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-15

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**