



Product Specification

Tube with false bottom CSF, 2.5 ml, (LxØ): 75 x 13 mm, round base, PP, Low Binding, cap enclosed, sterile, 1 piece(s)/blister

Page 1



Product description

Order number	63.614.625
Product description	Tube with false bottom CSF, nominal volume: 2.5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, Low Binding, with plastic label, label/print: transparent/blue, "CSF", fill mark, cap enclosed, blue, sterile, 1 piece(s)/blister, 100 piece(s)/carton

Product characteristics

Label/ Print	with plastic label
Colour of print/label	transparent/blue
Graduation	fill mark
Light protection	no
Application	collection and storage of diagnostic specimens
Example of use	Alzheimer's disease (AD) diagnostics
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	20

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-11-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Tube with false bottom CSF, 2.5 ml, (LxØ): 75 x 13 mm, round base, PP, Low Binding, cap enclosed, sterile, 1 piece(s)/blister

Page 2

Centrifugation time 10 Min
Shape of inner tube conical

Size

Sample volume 2.5 ml
Diameter 13 mm
Length of product 75 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap blue
Tube material Polypropylene (PP)

Purity & certification

Product category In-vitro diagnostic, CE
CE certified CE
Purity standard sterile
Sterilisation Gas sterilisation
Batched yes

Packaging

Minimum order qty. 100
Type of smallest blister
subpackaging
Piece(s) / inner box 10
Piece (s) / outer case 100
Piece(s) / pallet 12100
Depth of case 392 mm
Width of case 257 mm
Height of case 235 mm
Case volume 0.0237 cbm
Weight of product 0.0038 kg
Weight of case 0.82 kg
EAN of inner box 4038917896641

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-11-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Tube with false bottom CSF, 2.5 ml, (LxØ): 75 x 13 mm, round base, PP, Low Binding, cap enclosed, sterile, 1 piece(s)/blister

Page 3

EAN case

4038917896658

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-11-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT