



Product description

| | |
|---------------------|--|
| Order number | 65.1379.006 |
| Product description | Cap, blue, suitable for tubes Ø 12-17 mm, 1,000 piece(s)/bag |

Product characteristics

| | |
|--------------|------------------|
| Cap | standard cap |
| Suitable for | tubes Ø 12-17 mm |

Material & colours

| | |
|-------------------|----------------------------------|
| Product material | Low Density Polyethylene (LD-PE) |
| Colour of product | blue |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 10000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 10000 |
| Piece(s) / pallet | 240000 |
| Depth of case | 413 mm |
| Width of case | 257 mm |
| Height of case | 464 mm |
| Case volume | 0.0492 cbm |
| Weight of product | 0.0004 kg |
| Weight of case | 5.36 kg |
| EAN of inner box | 4038917128018 |
| EAN case | 4038917021548 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

