



## Product description

Order number	75.1356
Product description	Urine collection cup, collection and storage of urine, max. working volume: 500 ml, (ØxH): 80 x 148 mm, Ø opening: 80 mm, natural, material: PP, 250 piece(s)/case

## Product characteristics

Graduation	no
Application	collection and storage of urine
Volume of work	500 ml

## Size

Sample volume	500 ml
Max. volume	500 ml
Diameter	80 mm
Diameter of opening	80 mm
Height of product	148 mm
Length excluding cap	148 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine collection cup, 500 ml, (ØxH): 80 x 148 mm,**  
**PP, natural**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	natural

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

**Packaging**

Minimum order qty.	250
Type of smallest subpackaging	case
Piece(s) / inner box	250
Piece (s) / outer case	250
Piece(s) / pallet	5000
Depth of case	455 mm
Width of case	405 mm
Height of case	345 mm
Case volume	0.0636 cbm
Weight of product	0.03 kg
Weight of case	8.25 kg
EAN of inner box	4038917086318
EAN case	4038917072984

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**