




Product Specification
Urine container, 90 ml, (ØxH): 60 x 65 mm, PP,
with paper label, transparent



Product description

Order number	75.1435.004
Product description	Urine container, max. working volume: 90 ml, (ØxH): 60 x 65 mm, transparent, graduated, material: PP, screw cap, with paper label, sterile

Product characteristics

Label/ Print	with paper label
Graduation	yes
Closure type	screw cap
Volume of work	90 ml

Size

Sample volume	90 ml
Max. volume	90 ml
Diameter	60 mm
Height of product	65 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Urine container, 90 ml, (ØxH): 60 x 65 mm, PP,
with paper label, transparent

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	400
Piece(s) / inner box	400
Piece (s) / outer case	400
Piece(s) / pallet	7200
Case volume	0.1062 cbm
Weight of product	0.0201 kg
Weight of case	8.23 kg
EAN of inner box	4038917385220
EAN case	4038917385213

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

