



Product description

Order number	75.560.901
Product description	Urine container, max. working volume: 148 ml, (ØxH): 73 x 67 mm, white, material: PS

Product characteristics

Graduation	no
Volume of work	148 ml

Size

Sample volume	148 ml
Diameter	73 mm
Height of product	67 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 2500
Piece(s) / inner box 100
Piece(s) / case 2500
Piece(s) / pallet 90000
Depth of case 457 mm
Width of case 368 mm
Height of case 356 mm
Case volume 0.0599 cbm
Weight of product 0.0028 kg
Weight of case 8.16 kg
EAN of inner box 4038917318198
EAN case 4038917916547

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

