



### **Product description**

Order number	75.560.901
Product description	Urine container, max. working volume: 100 ml, (ØxH): 73 x 100 mm, white, material: PS

### **Product characteristics**

Graduation	no
Volume of work	100 ml

### **Size**

Sample volume	100 ml
Diameter	73 mm
Height of product	100 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Urine container, 100 ml, (ØxH): 73 x 100 mm, PS,**  
**white**

**Page 2**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

**Packaging**

Minimum order qty.                      2500  
Piece(s) / inner box                      100  
Piece(s) / case                              2500  
Piece(s) / pallet                            90000  
Case volume                                0,0762 cbm  
Weight of product                        0,0028 kg  
Weight of case                             8,16 kg  
EAN of inner box                         4038917318198  
EAN case                                    4038917318181

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-09-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**