



Product description

Order number	75.560.901
Product description	Urine container, max. working volume: 100 ml, (ØxH): 73 x 100 mm, white, material: PS

Product characteristics

Graduation	no
Volume of work	100 ml

Size

Sample volume	100 ml
Diameter	73 mm
Height of product	100 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Urine container, 100 ml, (ØxH): 73 x 100 mm, PS,
white**

Page 2

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 2500
Piece(s) / inner box 100
Piece(s) / case 2500
Piece(s) / pallet 90000
Case volume 0,0762 cbm
Weight of product 0,0028 kg
Weight of case 8,16 kg
EAN of inner box 4038917318198
EAN case 4038917318181

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT