



Product description

Order number	75.562.011
Product description	Container with screw cap, max. working volume: 100 ml, (ØxH): 57 x 76 mm, brown, with light protection, graduated, material: PP, screw cap, cap: brown, cap enclosed, 50 piece(s)/bag

Product characteristics

Graduation	yes
Type of cap	screw cap
Light protection	yes
Application	collection and storage of urine
Volume of work	100 ml

Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, with light protection, PP, brown

Page 2

Height of product 76 mm
Length excluding cap 72 mm

Material & colours

Product material Polypropylene (PP)
Colour of product brown
Cap material High Density Polyethylene (HD-PE)
Colour of cap brown

Purity & certification

Product category IVD product
CE certified CE
Batched yes

Packaging

Minimum order qty. 250
Type of inner packaging bag
Piece(s) / inner box 50
Piece(s) / case 250
Piece(s) / pallet 7500
Depth of case 443 mm
Width of case 323 mm
Height of case 330 mm
Case volume 0,0472 cbm
Weight of product 0,0153 kg
Weight of case 4,54 kg
EAN of inner box 4038917113618
EAN case 4038917014038

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT