



Product Specification
Container with screw cap, 100 ml, (ØxH): 57 x 76
mm, PP, transparent



Product description

Order number	75.562.300
Product description	Container with screw cap, max. working volume: 100 ml, (ØxH): 57 x 76 mm, transparent, graduated, material: PP, screw cap, cap: white, cap enclosed, 50 piece(s)/bag

Product characteristics

Graduation	yes
Type of cap	screw cap
Application	collection and storage of urine
Volume of work	100 ml

Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm
Height of product	76 mm
Length excluding cap	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Container with screw cap, 100 ml, (ØxH): 57 x 76
mm, PP, transparent

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	white

Purity & certification

Product category	IVD product
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	250
Type of inner packaging	bag
Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	7500
Depth of case	443 mm
Width of case	323 mm
Height of case	330 mm
Case volume	0,0472 cbm
Weight of product	0,0153 kg
Weight of case	4,59 kg
EAN of inner box	4038917113656
EAN case	4038917274739

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT