



Product description

Order number	75.570
Product description	Container, max. volume: 100 ml, (LxØ): 49 x 70 mm, Ø opening: 70 mm, white, cap: transparent, material: PS, snap-on lid, 50 piece(s)/tube bag

Product characteristics

Graduation	no
Closure type	snap-on lid
Volume of work	100 ml

Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	70 mm
Diameter of opening	70 mm
Height of product	49 mm
Length of product	49 mm
Length excluding cap	49 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Container, max. volume: 100 ml, (LxØ): 49 x 70 mm, PS, white

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	white
Colour of cap	transparent

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	tube bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	45000
Depth of box	320 mm
Width of box	75 mm
Height of box	75 mm
Depth of case	320 mm
Width of case	75 mm
Height of case	75 mm
Case volume	0.0317 cbm
Weight of product	0.0033 kg
Weight of case	3.71 kg
EAN of inner box	4038917082563
EAN case	4038917001649

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-08

This document was prepared by EDP support and is valid without signature.



SARSTEDT