



Product Specification
Medicine cup, 50 ml, (LxØ): 54 x 55 mm,
graduated, PP

Page 1



Product description

Order number	75.9924.522
Product description	Medicine cup, max. working volume: 50 ml, (LxØ): 54 x 55 mm, Ø opening: 55.4 mm, transparent, graduated, material: PP

Product characteristics

Graduation	yes
Application	Administration of medicines
Example of use	Oral administration of medications in liquid or tablet form
Volume of work	50 ml
Autoclavable (max. °C)	121 °C

Size

Sample volume	50 ml
Max. volume	50 ml
Diameter	55 mm
Diameter of opening	55,4 mm
Height of product	54 mm
Length of product	54 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-09-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Medicine cup, 50 ml, (LxØ): 54 x 55 mm,
graduated, PP

Page 2

Length excluding cap 55.4 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category Medical device
CE certified CE
Latex-free yes
Pyrogen-free no
Safety data sheet no

Packaging

Minimum order qty. 500
Piece(s) / inner box 500
Piece(s) / case 500
Piece(s) / pallet 80000
Depth of case 271 mm
Width of case 224 mm
Height of case 158 mm
Case volume 0.0099 cbm
Weight of case 3 kg
EAN of inner box 4038917352949
EAN case 4038917352932

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-09-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT