



### **Product description**

Order number	76.1340.500
Product description	Lid, for 30-ml medicine cup, material: PE, white, 750 piece(s)/bag

### **Product characteristics**

Example of use	Oral administration of medications in liquid or tablet form
Suitable for	30-ml medicine cup

### **Material & colours**

Product material	Polyethylene (PE)
Colour of product	white
Colour of cap	white

### **Purity & certification**

Product category	Medical device (MD)
CE certified	CE supplier
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## **Packaging**

Minimum order qty.	3750
Type of smallest subpackaging	bag
Piece(s) / inner box	3750
Piece (s) / outer case	3750
Piece(s) / pallet	180000
Depth of box	190 mm
Width of box	390 mm
Height of box	390 mm
Depth of case	395 mm
Width of case	190 mm
Height of case	400 mm
Case volume	0.03 cbm
Weight of product	0.0007 kg
Weight of case	2.97 kg
EAN of inner box	4038917127790
EAN case	4042836007149

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

