



Product description

Order number	76.571
Product description	Snap-on lid, for 100-ml container 75.570, material: PS, transparent, 100 piece(s)/tube bag

Product characteristics

Closure type	snap-on lid
Suitable for	100-ml container 75.570

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Snap-on lid, for 100-ml container 75.570, PS,
transparent

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	tube bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	32000
Depth of box	550 mm
Width of box	75 mm
Height of box	75 mm
Depth of case	550 mm
Width of case	75 mm
Height of case	75 mm
Case volume	0.0392 cbm
Weight of product	0.0009 kg
Weight of case	1.46 kg
EAN of inner box	4038917082570
EAN case	4038917001656

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

