



Product description

Order number	51.1703.009
Product description	Formalin Set, preparation: buffered formalin solution 450 ml, set components: Formalin system with 450 ml formaldehyde solution (4%) and 100 S-Monovettes 9 ml, 1 piece(s)/case

Product characteristics

Type of preparation	buffered formalin solution 450 ml
Label/ Print	with paper label
Design	S-Monovette® 9 ml
Application	Tissue fixation

Size

Width of product	98 mm
Length of product	330 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Formalin Set, buffered formalin solution 450 ml,
S-Monovette® 9 ml

Page 2

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	126
Depth of box	344 mm
Width of box	243 mm
Height of box	103 mm
Depth of case	350 mm
Width of case	250 mm
Height of case	115 mm
Case volume	0.0101 cbm
Weight of product	1.312 kg
Weight of case	1.77 kg
EAN of inner box	4038917336062
EAN case	4038917336055

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT