



Product description

Order number	51.1703.025
Product description	Formalin Set, preparation: buffered formalin solution 450 ml, set components: Formalin system with 450 ml formaldehyde solution (4%) and 50 S-Monovettes 25 ml, 1 piece(s)/case

Product characteristics

Type of preparation	buffered formalin solution 450 ml
Label/ Print	with paper label
Design	S-Monovette® 25 ml
Application	Tissue fixation

Size

Width of product	98 mm
Length of product	330 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	126
Depth of box	344 mm
Width of box	243 mm
Height of box	103 mm
Depth of case	355 mm
Width of case	255 mm
Height of case	123 mm
Case volume	0.0111 cbm
Weight of product	1.232 kg
Weight of case	1.69 kg
EAN of inner box	4038917336093
EAN case	4038917336086

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

