



## Product description

|                     |   |
|---------------------|---|
| Order number        | 51.1703.025   |
| Product description | Formalin Set, preparation: buffered formalin solution 450 ml, set components: Formalin system with 450 ml formaldehyde solution (4%) and 50 S-Monovettes 25 ml, 1 piece(s)/case |

## Product characteristics

|                     |                                   |
|---------------------|-----------------------------------|
| Type of preparation | buffered formalin solution 450 ml |
| Label/ Print        | with paper label                  |
| Design              | S-Monovette® 25 ml                |
| Application         | Tissue fixation                   |

## Size

|                   |        |
|-------------------|--------|
| Width of product  | 98 mm  |
| Length of product | 330 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



### **Purity & certification**

|                  |                                      |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD)            |
| CE certified     | CE - manufacturer's self-declaration |
| Batched          | yes                                  |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1             |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 1             |
| Piece(s) / pallet             | 126           |
| Depth of box                  | 344 mm        |
| Width of box                  | 243 mm        |
| Height of box                 | 103 mm        |
| Depth of case                 | 355 mm        |
| Width of case                 | 255 mm        |
| Height of case                | 123 mm        |
| Case volume                   | 0.0111 cbm    |
| Weight of product             | 1.232 kg      |
| Weight of case                | 1.69 kg       |
| EAN of inner box              | 4038917336093 |
| EAN case                      | 4038917336086 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

