



## Product description

Order number	51.1703
Product description	Formalin system, preparation: buffered formalin solution 450 ml, 1 piece(s)/inner box

## Product characteristics

Type of preparation	buffered formalin solution 450 ml
Label/ Print	with paper label
Application	Tissue fixation

## Size

Max. volume	450 ml
Width of product	98 mm
Length of product	330 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product in accordance with MDD or IVDD	IVD product
CE certified	CE
Lot or serial number	yes

### **Packaging**

Minimum order qty.	6
Type of inner packaging	inner box
Piece(s) / inner box	1
Piece(s) / case	6
Piece(s) / pallet	360
Depth of box	102 mm
Width of box	71 mm
Height of box	353 mm
Depth of case	363 mm
Width of case	325 mm
Height of case	164 mm
Case volume	0,0193 cbm
Weight of product	0,5302 kg
Weight of case	4,22 kg
EAN of inner box	4038917336031
EAN case	4038917336024

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

