



Product Specification

Gallipot, max. volume: 150 ml, (LxØ): 49 x 75 mm,
PP, translucent

Page 1



Product description

Order number	75.9924.517
Product description	Gallipot, max. volume: 150 ml, (LxØ): 49 x 75 mm, translucent, material: PP

Product characteristics

Application	Pathology storage container
Volume of work	150 ml
Autoclavable (max. °C)	121 °C

Size

Max. volume	150 ml
Diameter	75 mm
Width of product	75 mm
Height of product	49 mm
Length of product	49 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-01-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Gallipot, max. volume: 150 ml, (LxØ): 49 x 75 mm,
PP, translucent

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	translucent

Purity & certification

Latex-free	yes
Safety data sheet	no

Packaging

Minimum order qty.	1260
Piece(s) / inner box	1260
Piece (s) / outer case	1260
Piece(s) / pallet	30240
EAN of inner box	4038917155304
EAN case	4038917234917

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-01-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT