



Product description

| | |
|---------------------|---|
| Order number | 75.9924.518 |
| Product description | Gallipot, max. volume: 250 ml, (LxØ): 83 x 75 mm, translucent, material: PP |

Product characteristics

| | |
|------------------------|-----------------------------|
| Application | Pathology storage container |
| Volume of work | 250 ml |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|-------------------|---------|
| Max. volume | 250 ml |
| Diameter | 75 mm |
| Width of product | 77.2 mm |
| Height of product | 83 mm |
| Length of product | 83 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | translucent |

Purity & certification

| | |
|-------------------|-----|
| Latex-free | yes |
| Pyrogen-free | no |
| Safety data sheet | no |

Packaging

| | |
|------------------------|---------------|
| Minimum order qty. | 500 |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 16000 |
| Depth of case | 428 mm |
| Width of case | 372 mm |
| Height of case | 438 mm |
| Case volume | 0.0697 cbm |
| Weight of case | 8.7 kg |
| EAN of inner box | 4038917155311 |
| EAN case | 4038917039352 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

