



Product description

Order number	13.2406.601
Product description	V-Monovette® VD, for virus diagnostics, nominal volume: 4 ml, (LxØ): 75 x 13 mm, transparent, fill mark, push cap, cap light blue, with paper label, non-sterile, 100 piece(s)/bag

Product characteristics

Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	light blue
Graduation	fill mark
Cap	membrane cap
Type of cap	push cap
Light protection	no
Application	virus diagnostics
Base shape	round base
Volume of work	4 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
V-Monovette® VD, 4 ml, cap light blue, (LxØ): 75 x 13 mm, 100 piece(s)/bag

Size

Sample volume	4 ml
Diameter	13 mm
Length including cap	81 mm
Length excluding cap	75 mm

Material & colours

Product material	Polyethylene therephthalate (PET)
Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	light blue
Membrane material	rubber
Tube material	Polyethylene therephthalate (PET)

Purity & certification

Product category	IVD product
CE certified	CE
Purity standard	non-sterile
Safety data sheet	nein
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	39000
Weight of product	0,0056 kg
EAN of inner box	4038917790604
EAN case	4038917790611

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

