



Product description

Order number 65.590
Product description Screw cap, suitable for cold transport container

Product characteristics

Suitable for cold transport container

Material & colours

Colour of cap natural

Purity & certification

Product in accordance with MDD or IVDD no medical device | IVD
Lot or serial number no

Packaging

Minimum order qty. 2000
Piece(s) / inner box 100
Piece(s) / case 2000
Piece(s) / pallet 80000
Depth of case 428 mm
Width of case 388 mm
Height of case 225 mm
Case volume 0,0374 cbm
Weight of product 0,0024 kg
Weight of case 5,47 kg
EAN of inner box 4038917082723
EAN case 4038917001762

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

